














No.	Name	Picture	Model	Specification	Application
102	Magnetic Vibration & Heating Therapy Device		XY-K-CZR-III (LCD)	Combining the three treatment methods of magnetic therapy, vibration and heat therapy, can realize three treatments simultaneously; output channel: four channels; treatment mode: six treatment modes can be selected; Accessory: warm conductor, neck and shoulder conductor, joint conductor, separated conductor. Remarks: trolley is optional.	
103	Warm Magnetic Therapy Device		XY-K-GC-II	Thermal magnetic therapy for the spine, the temperature can be adjusted in five gears, standard treatment mat and matching treatment bed; treatment mat has over-temperature protection function, anti-scald; five treatment modes. Remarks: trolley is optional.	
	Accessory				
104	Pulse Magnetic Therapy Device	 	XY-K-MCC-II	Pulse current produces magnetic field with different frequency and waveform through the coil of electromagnet. Pulse magnetic therapy device has analgesic, anti-inflammatory effects, promote microcirculation, improve nutritional supply, reduce blood pressure and blood lipids.	
HEAT THERAPY DEVICE					
105	Bioceramics Hot Pack		XY-BC-II (1.75KG) (2.5KG)	Heat-resistant, when the temperature is 100°C, the bioceramics are not broken; press-resistant, the bioceramics can not break under the pressure of 10kg; the life-span of the bioceramics is 10 years.	
106	Bioceramics Hot Cushion		XY-BC-III (5kg, 10kg)	Heat-resistant, when the temperature is 100°C, the bioceramics are not broken; press-resistant, the bioceramics can not break under the pressure of 10kg; the life-span of the bioceramics is 10 years.	