










No.	Name	Picture	Model	Specification	Application
123	Gynecological Steam Therapy Device		XY-XZ-II	<p>Fumigation equipment specially designed for obstetrics and gynecology, outcropping full-body sitting fumigation cabin, microcomputer control operating system, double-layer glass fiber reinforced plastic material, easy to clean, and not easy to breed bacteria; designed on the principle of steam diversion, the airflow is soft, and the circulation rotation can treat all kinds of human body. Parts are fumigated, the effect is remarkable; it has automatic timing, automatic temperature control, automatic leakage protection, automatic anti-dry burning, automatic water replenishment, ozone disinfection, double temperature control protection, etc.</p> <p>Function: according to the ergonomic design, it can meet different heights, focus more on body curves, and avoid fatigue.</p>	
124	Steam Therapy Device		HYZ-IA	<p>Computer controlled, automatic temperature controlling, time setting, automatic leakage protection and dry burning prevention, LED display.</p>	
MOXIBUSTION THERAPY DEVICE					
125	Moxibustion Bed		XY-JLC-I	<p>Microcomputer controlled automatic ignition, easy to operate. Automatic translation: the moxibustion tray reciprocates during the treatment process to stimulate the acupoints more comprehensively. Automatic smoke exhaust: automatic smoke exhaust after treatment, which solves the problem of moxa smoke. Comfortable body position: supine governor meridian moxibustion, not easy to fatigue, enjoy the care of moxibustion in a relaxed way.</p>	
126	Moxibustion Bed		XY-JLC-III	<p>Microcomputer controlled automatic ignition, easy to operate. Automatic translation: the moxibustion tray reciprocates during the treatment process to stimulate the acupoints more comprehensively. Intelligent control: automatically adjust the temperature of the moxibustion tray according to the temperature inside the bed, so that the application part can be controlled at a constant temperature. Automatic smoke exhaust: automatic smoke exhaust after treatment, which solves the problem of moxa smoke. Comfortable body position: supine governor meridian moxibustion, not easy to fatigue, enjoy the care of moxibustion in a relaxed way.</p>	