

















No.	Name	Picture	Model	Specification	Application
HYDROTHERAPY DEVICE					
141	Ozone Hydrotherapy Bathtub		XY-SL-CIX	It is a fitness method, it uses the air and water of the most appropriate body as a therapeutic factor, foam board releases a certain pressure, temperature and gas in the water, and it naturally forms a large series of blisters and water droplets which generates many different level impacts on the human body skin, such as adsorption and massage. So in full of comfort and fun to achieve fitness, relieve fatigue of the muscles.	
142	Upper Limb Bathtub		XY-SL-BI	It is equipped with water jet massage function, through the spa mixed bubble massage on the upper extremities, LED colorful light, could adjust patient's psychological emotions. With the pump idle protection, to prevent water pump from damaged; electronic control disinfection system controls ozone generator to disinfect; electronic control drainage system is very convenient to drain away water; microcomputer control panel, more convenient to use and direct manipulation.	
143	Lower Limb Bathtub		XY-SL-BII	It is equipped with water jet massage function, massage could relax the lower leg and feet. LED colorful light, could adjust patient's psychological emotions, mental relaxation. With the pump idle protection, to prevent water pump from being damaged; electronic control disinfection system controls ozone generator to disinfect; electronic control drainage system is very convenient to drain away water; microcomputer control panel, more convenient to use and direct manipulation. Besides, temperature sensor monitors temperature change at any time.	
144	Lower Limb Bathtub		XY-SL-BIII	Six massage jets, massage lower leg with mix bubbles. LED colorful light, could adjust patient's psychological emotions, mental relaxation. With the pump idle protection, to prevent water pump being damaged; electronic control disinfection system controls ozone generator to disinfect; electronic control drainage system is very convenient to drain away water; microcomputer control panel, more convenient to use and direct manipulation. Besides, temperature sensor monitors temperature change at any time.	
145	Lower Limb and Spine Bathtub		XY-SL-BIV	Twelve massage jets, massage lower leg and spine with mix bubbles. LED colorful light, could adjust patient's psychological emotions, mental relaxation. With the pump idle protection, to increase pump service time; electronic control disinfection system controls ozone generator to disinfect; electronic control drainage system is very convenient to drain away water; microcomputer control panel is more convenient to use and direct manipulation. Besides, temperature sensor monitors temperature change at any time.	



No.	Name	Picture	Model	Specification	Application
146	Underwater Massage Device		XY-SL-BV	Equipped with pump idle protection system, increase service life. Pressure control valve adjusts water pressure to the appropriate range. Microcomputer control panel is more convenient to use and direct manipulation. Besides, temperature sensor monitors temperature change at any time. Equipped with wheel brakes and handles for easy movement and placement.	
147	Hydrotherapy and Electrotherapy Bathtub		XY-SL-BVII	<ul style="list-style-type: none"> • Liquid crystal display time and temperature, to simplify operation, besides own a variety of current mode. • There are 30 self-prescriptions in the device, operation more quick and convenient. • Intelligent heating mode meets the need of different treatment. • Temperature real-time monitoring system to ensure that the water temperature is comfortable; drainage function is controlled by button. 	
148	Butterfly Bathtub		XY-SL-CI	It provides sufficient motion space for the torso, arms and legs of the patient as special shape.	
149	Mobile Carrier		XY-SL-CF-1	<p>Butterfly shaped bathtub is suitable for patients of vascular accident, shoulder-hand syndrome, spinal cord injury, Parkinson's syndrome, fracture sequelae, post-burn rehabilitation.</p> <p>XY-SL-CF-1: used to move patients up or down horizontally. Space in bottom of the bath is sufficient for vertical movement of legs; and shallow both sides and underside of the head eliminate dead spots. And can stop at any height, easy to transfer patients from wheelchairs or other trolleys.</p>	
150	Mobile Carrier		XY-SL-CF-2	<p>XY-SL-CF-2: used to move patients up or down horizontally. Could stop at any height, easy to transfer patients from wheelchairs or other trolleys into butterfly bathtub.</p>	
151	Fixed Carrier		XY-SL-CF-3	<p>XY-SL-CF-3: used to move patients up or down horizontally. Could stop at any height, easy to transfer patients from wheelchairs or other trolleys into butterfly bathtub.</p>	









No.	Name	Picture	Model	Specification	Application
152	Full Body Massage Bathtub		XY-SL-CIII	<ul style="list-style-type: none"> • It is controlled by microcomputer and more convenient for operation; • Temperature real-time monitoring system to ensure that the water temperature is comfortable; • Electronic control disinfection system to ensure clean water; • Functions such as bubbles, vortex, constant temperature, surfing and ambient lights which can be chosen according to the clinical needs; • Adjustable water temperature in a range of 25°C-42 °C, and intelligent constant temperature; • Equipped with pump idling protection to extend the service life of the pump; • Warm and cold shower with convenient adjustment to provide a "pleasing" bath experience. 	
153	Full Body Automatic Massage Bathtub		XY-SL-CIV	<ul style="list-style-type: none"> • 256 water flow bubble nozzles, which can fully and partially open and close the water flow bubble function to provide the user with the maximum degree of freedom; • Equipped with electronic control disinfection system to control the ozone generator for disinfection; • Independent air pump, producing abundant air bubbles, modern automatic flushing and a drying system; • Equipped with pump idling protection to extend the service life of the pump; • Innovative electronic control disinfection, and a water descaling system to eliminate scale and disinfect water. 	
154	Rehabilitation Water Tank		XY-SL-CV	<ul style="list-style-type: none"> • Low-position entry, open-door bathtub and unique innovation to allow patient's easy access; • Hand-held sprinkler function: rotate the sprinkler to switch between the big flow mode and the massage mode; • Intelligent constant temperature function; • Three kinds of immersions: hydraulic massage, bubble massage, and soaking bath. 	
155	Rehabilitation Water Tank		XY-SL-CVI	<ul style="list-style-type: none"> • Hydraulic massage: multiple massage nozzle groups that allow choice of jet hydraulic massage and vortex hydraulic massage; • Thermostat device: intelligent thermostat device with automatic protection function; • Bubble massage: air is added to generate small warm bubbles, which move back and forth, so that oxygen-rich water fully covers in the body's soaked skin surface, making people feel happy and relaxed. 	



No.	Name	Picture	Model	Specification	Application
156	Full Body Electrotherapy Bath tub		XY-SL-CX	Set and display parameters through the touch screen to achieve simple and fast operation; with a variety of current modes to choose from, with multi-level adjustable intensity, for whole-body electric water bath; PTC static constant temperature system, built-in one-button quick drainage, equipped with high-pressure water massage device, can massage the parts of the human body, the massage pressure is adjustable, the massage pressure is displayed in real time, and it is equipped with a shower function, which is convenient for cleaning the bathtub. Equipped with leakage protection and overload protection to ensure the safe operation of the equipment.	
157	Ozone Hydrotherapy Bath tub		XY-SL-CIX	It adopts touch screen for operation and display, perfect combination of oxygen bath and bubble bath; humanized appearance design, automatic program design; constant temperature function; 7 fixed modes, 1 automatic cleaning mode and manual mode, and the massage force is adjustable; It has functions such as ozone sterilization. Equipped with leakage protection and overload protection to ensure the safe operation of the equipment. It can be used with the XY-SL-CF-1 mobile carrier to provide convenience for patients with limited mobility and reduce the work intensity of nursing staff.	
158	Smart Medicated Bath tub		XY-SL-CXI	Microcomputer control system, digital display of time and temperature, intelligent constant temperature mode to meet different treatment needs, temperature detection system to ensure constant water temperature. Equipped with leakage protection, overload protection and over-temperature protection to ensure the safe operation of the equipment. It can be used with the XY-SL-CF-1 mobile carrier to provide convenience for patients with limited mobility and reduce the work intensity of nursing staff.	
159	Dry Water Massage Table		XY-SL-CVIII	It uses buoyancy and softness and formlessness of water to lift up the body gently. Specially-made bed surface can not only bear the pressure of all patients' weight, but also can effectively relay water massage force. The bed surface can fit the curve of human body, so that the cervical spine, lumbar spine and other parts will do not suspend in air, and the body parts come under uniform stress, the patient can feel that his body is relaxed naturally. Then the force is acted upon human body for massage therapy by means of the water flow force, the motion mode and motion position of the microcomputer-controlled high-energy nozzle.	
			XY-SL-CVIII (deluxe)		



No.	Name	Picture	Model	Specification	Application
160	Dry Water Massage Table		XY-SL-CXIII	High-pressure water impacts the bottom of the silicone bed surface of the massage table for body massage, soothes muscles, relieves pressure, eliminates pain, smooths veins and lymphatic vessels, relaxes subcutaneous tissue, promotes blood circulation, and stimulates metabolism. It adopts a movable flat panel for operation and display. It has 13 massage modes, 2 independent shoulder nozzles can be turned on or off independently, 6 kinds of massage parts to choose from, 10 levels of pressure adjustment, 10 levels of speed adjustment, temperature setting and display, time setting and display, water-cooled heat dissipation, automatic exhaust, height setting and other functions, the equipment can perform massage with up to 6 nozzles at the same time. At the same time, it has multiple functions such as music, aroma, air purification, ambient light, and foot massage. It has functions such as leakage protection, over-temperature protection, overload protection, and dry-burning protection to ensure the safe operation of the equipment.	
161	Aquatic Underwater Treadmill		XY-SL-CXII (without filter)	This system has the reverse vortex water bubble, wide adjustable motor, screen visual feedback, suction cup type handrail, rapid water circulation disinfection device and other practical functions to make sure each patient has proper treatment and better training effect.	
162	Lifting Chair		XY-SL-CF-5	The mobile lift chair adopts the electric lifting column structure to realize the height lifting of the patient.	
163	Underwater Treadmill		XY-SL-CF-6	Dimension: 1230mm*620mm*1540mm; composed of hand lever, support frame, handle, rolling shaft, etc; 19 rolling axes; an adsorption device, with anti-slip function.	
164	Underwater Static Bike		XY-SL-CF-7	Dimension: 965mm*551mm*1485mm; composed of hand lever, support frame, transmission mechanism, etc.	
165	Underwater Bike		XY-SL-CF-9	With the buoyancy of water, resistance, and walking training with the strength of lower limbs, it is suitable for endurance training, gait training, and training of lower limb joint activities.	