




















No.	Name	Picture	Model	Specification	Application
<b>PHOTOTHERAPY DEVICE</b>					
15	High Power Infrared Therapy Device		XY-K-GNHW-II	It uses a halogen light source as the luminous body, after the light is filtered by a special liquid filtering system and a filter, the light that is prone to thermal effects on the skin is greatly filtered out, and only the light that the patient can tolerate and has great therapeutic value is retained. A high-energy light wave (wavelength range 560nm~1400nm) that can easily penetrate 10cm under the human skin to form a deep treatment.	
16	Medical Infrared Thermal Imager		XY-RX-I	It is used to check acute and chronic pain in bone joints, meridian, muscle and soft tissue in whole body.	
17	Infrared Therapy Device		XYG-500IB	With large radiation range, stable performance, and outstanding treatment effects. Beautiful and convenient trolley design, green treatment, no infection and no complications.	
18	Infrared and Laser Therapy Device		XYG-500IIB	With the design of punctiform radiator and planar radiator, LCD display, with an emergency switch.	
19	Infrared and Laser Therapy Device		XYG-500IIIB	With the design of punctiform and planar radiator, LCD display, with an emergency switch, multi-joint mechanical arm and key switch.	
20	Infrared and Laser Therapy Device		XYG-500IVB	With the design of punctiform and planar radiator; stepwise treatment; triple area and punctiform radiator.	
21	Ultraviolet Radiation Therapy Device		XY-K-ZWX-II	With the wavelength of 254 nm, ultraviolet ray has good drying, sterilization, anti-inflammatory effects. It can accelerate blood circulation, analgesia, promote epithelial tissue regeneration, without adverse response, simple and safe treatment.	



No.	Name	Picture	Model	Specification	Application
22	UV Sterilizer		XDD-I	<ol style="list-style-type: none"> <li>1. Dual-channel ultraviolet lamp output, high-purity quartz lamp, good sterilization effect.</li> <li>2. Height-adjustable bracket, trolley design, easy to move.</li> <li>3. The chassis of the trolley adopts X-shaped sheet metal design, which is light in weight and elegant in appearance.</li> <li>4. Easy to install and use, takes up little space.</li> </ol>	
23	UV Sterilizerr		XDD-II	<ol style="list-style-type: none"> <li>1. The base of the sterilizer can rotate flexibly, and the height can be automatically retracted.</li> <li>2. Exquisite trolley design, easy to move.</li> <li>3. The UV switch can be controlled by a wireless remote control to prevent harm to the human body.</li> <li>4. It can automatically turn and stop at a fixed time.</li> <li>5. Human body infrared sensor safety protection will automatically shut down when it detects that someone enters the disinfection area.</li> <li>6. Three-way ultraviolet lamp output, high-purity quartz lamp, better sterilization effect.</li> </ol>	
24	UV Sterilizer		XDD-III	<ol style="list-style-type: none"> <li>1. 222nm and 254nm integrated design, to achieve a safe environment for human-computer coexistence.</li> <li>2. 360° full-angle human body induction linkage control.</li> <li>3. Portable regular hexagon shape design, more stable and reliable, easy to move.</li> <li>4. With remote control and operation voice broadcast, alarm, the operation is clear and distinct.</li> <li>5. Equipped with indicator lights of different colors.</li> <li>6. The cumulative working hours statistics function makes maintenance more convenient.</li> </ol>	
25	Infrared Thermometer		XY-ETW-I	<ol style="list-style-type: none"> <li>1. Two modes: surface temperature mode and body temperature mode;</li> <li>2. Memory query: 32 sets of measurement data can be memorized for easy query;</li> <li>3. temperature unit can be selected from °C or °F;</li> <li>4. Automatic shutdown.</li> </ol>	
26	Infrared Polarized Light Therapy Device		XY-K-PZG-II	<p>Infrared polarized light produces intense light and acupuncture moxibustion effects, which has an effect on the body's nervous system, circulatory system, cardiovascular, digestive system, endocrine system and immune system adjustment, which changes the pathophysiologic process of the body, restore physiological balance and maintain a stable internal environment.</p>	
27	Infrared Polarized Light Therapy Device		XY-K-PZG-III		