






No.	Name	Picture	Model	Specification	Application
<b>PARAFFIN WAX THERAPY DEVICE</b>					
107	Paraffin Wax Therapy Device		XYL-I	Anti-electric leakage function and leakage protection device; nohydration wax technology; automatic temperature control device, keep the temperature at 58°C-60°C; melt wax: 10kg.	
108	Paraffin Wax Therapy Device		XYL-II	Anti-electric leakage function and leakage protection device; nohydration wax technology; automatic temperature control device, keep the temperature at 58°C-60°C; melt wax: 30 kg.	
109	Paraffin Wax Therapy Device		XYL-V	Anti-electric leakage function and leakage protection device; nohydration wax technology; automatic temperature control device, keep the temperature at 58°C-60°C; melt wax: 30kg.	
110	Paraffin Wax Therapy Device		XYL-VIIA XYL-VIIC	The newly upgraded intelligent control, advanced and very user-friendly functions such as voice broadcast, power-off memory, intelligent ultraviolet disinfection, and automatic switching of wax making procedures;  Three cake making modes: Normal wax making, rapid wax making, and reservation wax making to meet the needs of different situations; 5 electrical safety protection functions, overheating protection, current overload protection, fault self-checking alarm and other functions; ergonomic design, frosted stainless steel oversized worktable; with various wax plates, to make various sizes of wax cakes. Wax cake can be up to 12-36 layers, and a maximum of 28-84 wax cakes can be produced.	
111	Paraffin Wax Therapy Device		XYL-VIID XYL-VIIF		
112	Operating Table		XY-LF-I		
<b>STEAM THERAPY DEVICE</b>					
113	Steam Therapy Device		HYZ-IIA	Microcomputer control operating system, digital display window; with timing, automatic temperature control, automatic leakage protection, automatic anti-dry, ozone disinfection function; foot massage and Chinese medicine liquid atomization function, which can be used for whole body fumigation and reproductive system fumigation.	



No.	Name	Picture	Model	Specification	Application
114	Steam Therapy Device		HYZ-IIIB (touch screen display)	The herbs steam temperature and time are controlled by the microcomputer, with overheat protection function, the cabin is made of glass fiber reinforced plastics, the LCD screen can display temperature, time, water level, the intelligent software can set the fumigation programs, dynamic monitoring temperature, time and water level, with leakage protection function, automatic temperature controlling function, automatic setting time function, self-sterilization to avoid the cross infection.	
115	Steam Therapy Device		HYZ-IIIE	The height of the seat is adjustable, it is suitable for the persons who are 1.5-2.0m tall, the person can sit, almost lie down and lie down, with oxygen supply, the patient can breathe the oxygen if their chest does not feel well, with atomization function, separated console, electromechanical separation make the machine safe and reliable, with MP3 player and multifunctional shower head.	
116	Steam Therapy Device		HYZ-IIH	Glass fiber reinforced plastics main body, and can choose multiple fumigation method; three heating zones can be used separately to make the treatment having more benefits.	
117	Steam Therapy Device		HYZ-IIIC (with trolley)	Microcomputer control, can set temperature automatically, can imitate acupuncture and massage, impacting treatment, independent console, with atomization function(for the machine with console).	
118	Steam Therapy Device		HYZ-IIID (for foot and hand)	This machine is for foot and hand fumigation, three persons can do the fumigation treatment at the same time, microcomputer control, automatic time controlling and temperature controlling, automatic burning prevention, with ozone sterilization function and MP3 player.	
119	Steam Therapy Device		HYZ-IE	Brand new upgrade, large-screen LCD touch control, intelligent skin temperature control, multi-channel heat dissipation system, three-dimensional and multi-angle nozzles, frequency conversion control heating, and multiple safety protection measures against over-temperature or over-pressure or leakage, liquid shortage and dry burning. For effective protection, it has many practical functions such as pre-heating decoction, pre-heating insulation, intelligent discharge of waste liquid, and it can monitor skin temperature (patent) by infrared in real time to prevent patients from being scalded during treatment.	



No.	Name	Picture	Model	Specification	Application
120	Steam Therapy Device		HYZ-IB	Brand new human-computer interactive system. Dissolve TMC by microcomputer. Large-screen LCD touch control. Several safety protection measures against over-temperature, over-pressure, leakage, liquid shortage and dry burning. For effective protection, it has many practical functions such as pre-heating decoction, pre-heating insulation, intelligent discharge of waste liquid, and it can monitor skin temperature (patent) by infrared in real time to prevent patients from being scalded during treatment.	
121	Steam Therapy Device		HYZ-IIID	Large-space luxury outcropping full-body sitting fumigation cabin, microcomputer-controlled operating system, large-screen capacitive touch; double-layer glass steel material, easy to clean, not easy to breed bacteria; can be used for outcrop full-body fumigation; innovative multi-point hydro massage on the back, fast eliminate fatigue and relax tension; the design of the diversion principle, the airflow is soft, and the circular rotation can fumigate the joints of the human body, with remarkable effects; it has automatic timing, automatic temperature control, automatic leakage protection, automatic anti-dry heating, automatic water intake, fog functions such as disinfection, ozone disinfection, hydromassage, MP3 playback, dual temperature control and protection; according to the ergonomic design, the height of the seat can be adjusted to meet different heights, and more focus points can fit the body curve to avoid fatigue; it is equipped with a multi-functional shower you can bathe as you like.	
122	Steam Therapy Device		HYZ-IIY	Newly upgraded, full-featured AI chip, voice prompts, microcomputer touch screen display control operating system; intelligent design, patented technology, one-key completion of ozone flushing, traditional Chinese medicine fumigation, water flushing, drying and other functions, easy to operate; low power consumption, fast heating, multi-mode operation, automatic, manual, and individual mode treatment to meet the needs of different patients; functions such as preheating and heat preservation, automatic water intake, anti-dry heating, ozone disinfection, bluetooth music, seat heating, dual temperature protection, etc.; high-precision backwash filtration ensures the smooth operation of the equipment, and the one-key self-cleaning function reduces the tedious cleaning steps of medical staff; the independent control console better protects the privacy of patients.	