



No.	Name	Picture	Model	Specification	Application
127	Moxibustion Bed		XY-JLC-IV	LCD display, tablet APP control, automatic ignition, easy to operate. Automatic translation: the moxibustion tray reciprocates during the treatment process to stimulate the acupoints more comprehensively. Intelligent control: automatically adjust the height of the moxibustion tray according to the temperature inside the bed, so that the applied part can be controlled at a constant temperature. Automatic smoke exhaust: Automatic smoke exhaust after treatment, which solves the problem of moxa smoke. Comfortable body position: supine governor meridian moxibustion, not easy to fatigue, enjoy the care of moxibustion in a relaxed way. Cabin moxibustion, better heat preservation, more concentrated treatment effect.	
128	Moxibustion Bed		XY-AJC-I	16 moxibustion boxes can be selected according to the needs of patients to choose parts, accurate treatment; equipped with storage shelves under the bed, easy to use; two sections bed board, humanized design, adjustable tilt angle.	
129	Moxibustion Bed		XY-AJC-III	Microcomputer control auto-ignition, easy to operate; intelligent regulation: ensure the efficacy; automatic smoke exhaust, automatic smoke exhaust after treatment.	
HYDROCOLLATOR					
130	Hydrocollator		XY-SRF-I	The device is made of stainless steel, thermostatical control, easy to operate, the knob to adjust the temperature, 13 pieces of hot packs, 10 standard hot packs, 1 lengthen hot pack, 1 hot pack for neck, 1 hot pack for hand.	
131	Hydrocollator		RH-SR-II	<ol style="list-style-type: none"> 1. Constant temperature: 0~99°C adjustable; 2. Over-temperature protection and low water level, alarm to prevent drying heating; 3. Material: 304 stainless steel and fiberglass; 4. Knob controller ; 5. 38L; 6. Stationary type; Eight standard packs.	